



Implant Guide Kit

- Used for checking crown size and angulation of implant placement
- Provides placement guidance with one drilling procedure
- Guides angulation for 4 surfaces (mesial, distal, lingual and buccal)

Direction for use





COMPONENTS







Screw Remover Kit

- Removes broken screws to save the implant
- · Works for all major implant brands











COMPONENTS







Fixture Remover Kit

- Removes fixture to preserve the surgery area to place a new one
- Fixture can be removed by the Implant Forcep or the Fixture Remover
- In case of damage in the fixture connection, Trephine Drills are available to remove the fixture

Direction for use









Trephine CASE



COMPONENTS








Bone Expander Kit

- Expands narrow bone ridge or bad bony tissue
- Reversed thread system on the Expander provides better stability of the fixture
- Helps initial fixation on low density bones by compressing the bone outward

Direction for use











Ratchet

Crest Trimmer

Sawing

Expanding

Insert Implant

COMPONENTS



Universal Prosthetic Kit



- Consists of drivers that are compatible with all major implant brands
- Reduces carrying many different implant systems
- Helps tighten clinical and prosthetic screws to get the correct torque value

COMPONENTS

PRODUCT	COMPATIBLE IMPLANT	PICTURE	SIZE
1.2 HEX DRIVER	SIS / Neobiotech / Biomet 3i / Dentsply - Friadent, Xive / Osstem, DIO / Megagen		
1.25 HEX DRIVER	Biohorizons / Zimmer / Astra / Implant Direct Legacy / Camlog / MIS / Dentis / Dentiu		10415
SCS DRIVER	Straumann / Implant Direct Replant		10/15mm
UNIGRIP DRIVER	Nobel Biocare, Implant Direct Swishplait / Southern Implants / Neoss		
RATCHET			

